

# Cool Beauty Instrument



User Manual

**SK-0814**

Please read this manual carefully before operation

## Focus on quality Revolve about creation

Mystery for tightening skin

——Cool beauty instrument

### Cool beauty instrument

We devote ourselves to the development of beauty products for ladies, and are striving to create a sophisticated product cool beauty instrument in the beauty filed with professional quality, dedicated spirit and advanced technologies, so as to help you get rid of the age trail and renew your bright watery skin!

Cool beauty instrument is a new beauty product developed and produced by our company, authorized by Skinay Group (France) Co., Ltd..

We provide with a unique professional beauty production development team. After diligent research, cool beauty instrument is now proudly created. It is smart and easy to operate. Its first class quality and comfortable feelings make it more popular among people all over the world seeking for beauty, which makes us believe Skinary will change the beauty world!

## Cool beauty instrument

Do you still worry about slack skin?  
Are you still being bothered by rough pore?  
Don't worry! Cool beauty instrument will be helpful to you!

Cool beauty instrument is a new unique equipment by our company for tightening skin. It can quickly minimize the pores, eliminate the wrinkles and tighten the skin. It is smart and easy to use, and it is also convenient for daily skin care. It will make your skin better and brighter after long-time use.



## Do you know

Along with the lapse of time, the original collagen of skin will lose gradually. This will make the skin slacker, and the pores may be coarser. Hence, it is essential for ladies to tighten the skin and minimize the pores.



## Do you know

Expert suggestions: The optimum method of minimizing the pores and tightening the skin is to make ice care. According to the principle of thermal expansion and contraction, ice care can not only minimize the pores immediately but also decrease the surface temperature of skin quickly. In this way, we can effectively control the oily skin, and thus reach the best tightening effect.



## Freeze the skin quickly in 30s

Cool beauty instrument can cool down to 8°C in 30s, and can tighten the skin and minimize the pore. At the same time, it can also stimulate the collagen in skin, so as to increase elasticity and brightness of the skin!



Ice skin  
treatment →




Tightening

Smoothing

Elasticity

Brightness

## Product structure



## Product trait

### 1. Structure design

Cool beauty instrument is smart, simple and convenient to carry.


### 2. Function trait

It can cool down to 8℃ quickly and is applicable to each part of the face.

### 3. Care effect


The unique triangle design is specially designed for tightening the two side of nasal ala and T part, so as to care the skin completely.

## Product function



## Operation procedures

- 1.Remove the protective cover of ice skin treatment surface.
- 2.Power on.
- 3.Open the switch.
- 4.Take 90s to completely cool down.



## Use frequency

It can not be overused

It is suggested to use twice to three times each week at beginning, then adjust the use frequency and time according to the conditions of your skin (once per day at most).

## Use method

- Encircle around the face at the speed of 30s per circle, including the cheek and frontal angle.
- Pay more attention to some parts, such as the nasal ala.

## Care procedures

- 1.Move slowly in mode of circling.
- 2.Move from the internal to external and from lower part to upper part.
- 3.Symmetrical at both the left and right sides

## Tips

- Clean face with warm water first, then treat the skin with cool beauty instrument.
- Use astringent or emulsion on the face after completion of treatment.
- It is feasible to retreat according to the conditions of your skin.
- It is suggested to use twice to three times each week.

## Technical parameters and production specifications

Product name: Cool beauty instrument

Product model: SK-0814

Rated voltage: 100-240V

Rated power: 3.5V/3.5A

Product Dimensions: 18×4.5×3.5cm

Product Weight: 137g

### Notes:

The following phenomena are normal

1. The temperature of treatment surface will rise quickly after use and power dump;
2. There is light noise generated during application.

## Troubleshooting

Symptom	Cause	Solution
Ice skin beauty instrument can not cool down.	The power of main frame is not opened.	Start the power of main frame.
	The air inlet and outlet are blocked.	Hold the main frame correctly, take care not to block the air inlet and outlet.
	The power plug is not correctly inserted.	Re-insert the power plug.
	The power cable is cut off.	Check whether the power cables are normal or not.

### Notes:

The following phenomena are normal

1. The temperature of treatment surface will rise quickly after use and power dump;
2. There is light noise generated during application.

## Precautions

- 1、 This product shall be used under ambient temperature of 15 - 40°C .
- 2、 Check whether the frame and plug are damaged before use. If it is true, do not use.
- 3、 Keep the hands dry during energization.
- 4、 Do not treat the area around the eye during use.
- 5、 If there are erythema, hyperemia or other abnormalities, stop using immediately.
- 6、 Do not treat the injured skin with this apparatus.
- 7、 Do not use in damp conditions.
- 8、 Do not place this product at the places where the children can touch it easily.

## Cleaning and maintenance

- 1.Close the switch and power off firstly after use each time.
- 2.Swab the frame clean with wet towel after use.
- 3.Do not wash the whole product in water.
- 4.Store it in dry conditions after cleaning.

